

Provox LaryTube - When and How to Use

The Provox LaryTube is a cannula designed to hold Provox HMEs and Provox FreeHands HMEs and is placed in the stoma to prevent it from shrinking. It can also be used instead of an adhesive when the skin is irritated and it can be used during radiotherapy on the skin around the stoma when it is not possible to use an adhesive.


The Provox LaryTube is available in three versions:

- Fenestrated version – For voice prosthesis users. Designed with holes so that air can go through the voice prosthesis. Can be attached with a Provox TubeHolder or Provox LaryClips.
- Standard version – made for use with or without a voice prosthesis. Custom punctures can be made using the Provox Fenestration Punch.
- Blue Ring version – Can be attached with a Provox Adhesive. Recommended when the LaryTube is only used during the night. Punctures can be made with the fenestration punch for voice prosthesis users.



1. Sizing of the Provox LaryTube

- The Provox LaryTube comes in 4 widths. The tracheostoma should be just slightly bigger than the Provox LaryTube.
- The Provox LaryTube comes in 3 lengths.
- It is important that the Provox LaryTube you use is of the correct type, width and length, and fits well into your stoma without causing irritation. Your clinician can help you with determining which LaryTube is the best solution for you.



TIP: Be careful when inserting the Provox LaryTube that the end of the tube does not “hook” onto the prosthesis and cause it to move. Contact your clinician if you experience irritation from the tube, or if you feel the tube has become too large or too small. You may need a different size.

2. Inserting the Provox LaryTube


- Lightly moisten the tube with water soluble lubricant.
- Insert the Provox LaryTube into the tracheostoma.
- Attach the Provox LaryTube with Provox LaryClips, a Provox TubeHolder, or a Provox Adhesive.
- Insert Provox HME or FreeHands HME into the cassette holder.



Provox LaryTube with
Provox LaryClips

3. Cleaning the LaryTube

- Clean the LaryTube using a Provox TubeBrush. The Provox Brush can be used to clean the punctures of a fenestrated LaryTube.
- Disinfect the LaryTube in accordance with the LaryTube manual.



Provox Brush



Provox Tube Brush